# TruSight™ Oncology 500 ctDNA v2

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製品名	カタログ番号
ローカルのバリエントレポート作成	
Illumina DRAGEN Server v4	20051343
Illumina DRAGEN Server Installation	20031995
Illumina DRAGEN Server v4 Support Plan	20085832
Field Delivered Applications Training	15032919
クラウドベースのバリエントレポート作成	
ICA Basic Annual Subscription	20044874
ICA Professional Annual Subscription	20044876
ICA Enterprise Annual Subscription	20038994
ICA Enterprise Compliance Add-on	20066830
Subscription ICA Training and Onboarding	20049422
バリエント解釈	
Illumina Connected Insights – Annual Subscription	20090137
Illumina Connected Insights-Research – Annual Subscription	20112516
Illumina Connected Insights – Oncology Genome Equivalent Samples (VCF)	20090138
Illumina Connected Insights Training – Remote	20092376
Informatics Professional Services	20071787
クラウドストレージ	
Illumina Analytics – 1 iCredit	20042038
Illumina Analytics Starter Package – 1,000 iCredits	20042039
Illumina Analytics – 5,000 iCredits	20042040
Illumina Analytics – 50,000 iCredits	20042041
Illumina Analytics – 100,000 iCredits	20042042

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1. Lim C, Tsao MS, Le LW, et al. Biomarker testing and time to treatment decision in patients with advanced nonsmall-cell lung cancer. *Ann Oncol*. 2015;26(7):1415-1421. doi:10.1093/annonc/mdv208
2. Yu TM, Morrison C, Gold EJ, Tradonsky A, Layton AJ. Multiple Biomarker Testing Tissue Consumption and Completion Rates With Single-gene Tests and Investigational Use of OncoPrint Dx Target Test for Advanced Non-Small-cell Lung Cancer: A Single-center Analysis. *Clin Lung Cancer*. 2019;20(1):20-29.e8. doi:10.1016/j.clcc.2018.08.010
3. Reitsma M, Fox J, Borre PV, et al. Effect of a Collaboration Between a Health Plan, Oncology Practice, and Comprehensive Genomic Profiling Company from the Payer Perspective. *J Manag Care Spec Pharm*. 2019;25(5):601-611. doi:10.18553/jmcp.2019.18309
4. Zehir A, Benayed R, Shah RH, et al. Mutational landscape of metastatic cancer revealed from prospective clinical sequencing of 10,000 patients [published correction appears in Nat Med. 2017 Aug 4;23(8):1004]. *Nat Med*. 2017;23(6):703-713. doi:10.1038/nm.4333
5. Kopetz S, Mills Shaw KR, Lee JJ, et al. Use of a Targeted Exome Next-Generation Sequencing Panel Offers Therapeutic Opportunity and Clinical Benefit in a Subset of Patients With Advanced Cancers. *JCO Precis Oncol*. 2019;3:PO.18.00213. Published 2019 Mar 8. doi:10.1200/PO.18.00213
6. Drilon A, Wang L, Arcila ME, et al. Broad, Hybrid Capture-Based Next-Generation Sequencing Identifies Actionable Genomic Alterations in Lung Adenocarcinomas Otherwise Negative for Such Alterations by Other Genomic Testing Approaches. *Clin Lung Cancer*. 2015;21(16):3631-3639. doi:10.1158/1078-0432.CCR-14-2683
7. Buchhalter I, Rempel E, Endris V, et al. Size matters: Dissecting key parameters for panel-based tumor mutational burden analysis. *Int J Cancer*. 2019;144(4):848-858. doi:10.1002/ijc.31878
8. Chalmers ZR, Connelly CF, Fabrizio D, et al. Analysis of 100,000 human cancer genomes reveals the landscape of tumor mutational burden. *Genome Med*. 2017;9(1):34. Published 2019 Mar 8. doi:10.1186/s13073-017-0424-2
9. Hagemann IS, Devarakonda S, Lockwood CM, et al. Clinical next-generation sequencing in patients with non-small cell lung cancer. *Cancer*. 2015;121(4):631-639. doi:10.1002/cncr.29089
10. Bettgowda C, Sausen M, Leary RJ, et al. Detection of circulating tumor DNA in early- and late-stage human malignancies. *Sci Transl Med*. 2014;6(224):224ra24. doi:10.1126/scitranslmed.3007094
11. Saarenheimo J, Eigelienė N, Andersen H, Tirola M, Jekunen A. The Value of Liquid Biopsies for Guiding Therapy Decisions in Non-small Cell Lung Cancer. *Front Oncol*. 2019;9:129. doi:10.3389/fonc.2019.00129
12. Rapisuwon S, Vietsch EE, Wellstein A. Circulating biomarkers to monitor cancer progression and treatment. *Comput Struct Biotechnol J*. 2016;14:211-222. Published 2016 Jun 1. doi:10.1016/j.csbj.2016.05.004
13. Aggarwal C, Thompson JC, Black TA, et al. Clinical Implications of Plasma-Based Genotyping With the Delivery of Personalized Therapy in Metastatic Non-Small Cell Lung Cancer. *JAMA Oncol*. 2019;5(2):173-180. doi:10.1001/jamaoncol.2018.4305
14. Leigh NB, Page RD, Raymond VM, et al. Clinical Utility of Comprehensive Cell-free DNA Analysis to Identify Genomic Biomarkers in Patients with Newly Diagnosed Metastatic Non-small Cell Lung Cancer. *Clin Lung Cancer*. 2019;25(15):4691-4700. doi:10.1158/1078-0432.CCR-19-0624
15. Palmero R, Taus A, Viteri S, et al. Biomarker Discovery and Outcomes for Comprehensive Cell-Free Circulating Tumor DNA Versus Standard-of-Care Tissue Testing in Advanced Non-Small Cell Lung Cancer. *JCO Precision Oncology*. 2021;5:93-102. doi:10.1200/PO.20.00241
16. Mack PC, Banks KC, Espenschied CR, et al. Spectrum of driver mutations and clinical impact of circulating tumor DNA analysis in non-small cell lung cancer: Analysis of over 8000 cases. *Cancer*. 2020;126(14):3219-3228. doi:10.1002/cncr.32876
17. Rolfo C, Mack P, Scagliotti GV, et al. Liquid Biopsy for Advanced NSCLC: A Consensus Statement From the International Association for the Study of Lung Cancer. *J Thorac Oncol*. 2021;16(10):1647-1662. doi:10.1016/j.jtho.2021.06.017
18. Mosele F, Remon J, Mateo J, et al. Recommendations for the use of next-generation sequencing (NGS) for patients with metastatic cancers: a report from the ESMO Precision Medicine Working Group. *Ann Oncol*. 2020;31(11):1491-1505. doi:10.1016/j.annonc.2020.07.014
19. Pascual J, Attard G, Bidard FC, et al. ESMO recommendations on the use of circulating tumour DNA assays for patients with cancer: a report from the ESMO Precision Medicine Working Group. *Ann Oncol*. 2022;33(8):750-768. doi:10.1016/j.annonc.2022.05.520
20. Illumina. TruSight Oncology 500 ctDNA data sheet. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/trusight-oncology-500-ctdna-data-sheet-m-gl-00843/trusight-oncology-500-ctdna-data-sheet-m-gl-00843.pdf](https://www.illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/trusight-oncology-500-ctdna-data-sheet-m-gl-00843/trusight-oncology-500-ctdna-data-sheet-m-gl-00843.pdf). Accessed September 20, 2023.

21. Roche. AVENIO ctDNA Analysis Kits. Sequencing. [sequencing.roche.com/us/en/products/product-category/avenio-ctdna-analysis-kits.html](https://sequencing.roche.com/us/en/products/product-category/avenio-ctdna-analysis-kits.html). Accessed September 20, 2023.
22. Tempus. Tempus xF. [www.tempus.com/oncology/genomic-profiling/xf/](http://www.tempus.com/oncology/genomic-profiling/xf/). Accessed September 20, 2023.
23. Personal Genome Diagnostics. PGDx elio plasma complete. <https://www.personalgenome.com/products/elio-plasma-complete>. Accessed September 20, 2023.
24. Illumina. TruSight Oncology UMI Reagents data sheet. [illumina.com/content/dam/illumina-marketing/documents/products/datasheets/trusight-oncology-umi-reagents-datasheet-1000000050425.pdf](https://illumina.com/content/dam/illumina-marketing/documents/products/datasheets/trusight-oncology-umi-reagents-datasheet-1000000050425.pdf). Accessed September 20, 2023.
25. Illumina. Sequencing accuracy with Unique Molecular Identifiers. [illumina.com/techniques/sequencing/ngs-library-prep/multiplexing/unique-molecular-identifiers.html](https://illumina.com/techniques/sequencing/ngs-library-prep/multiplexing/unique-molecular-identifiers.html). Accessed September 20, 2023.
26. Socea JN, Stone VN, Qian X, Gibbs PL, Levinson KJ. Implementing laboratory automation for next-generation sequencing: benefits and challenges for library preparation. *Front Public Health*. 2023;11:1195581. doi:10.3389/fpubh.2023.1195581
27. Data on file. Illumina, Inc. 2024.
28. Beroukhim R, Mermel CH, Porter D, et al. The landscape of somatic copy-number alteration across human cancers. *Nature*. 2010;463(7283):899-905. doi:10.1038/nature08822
29. Illumina. Analysis of TMB and MSI status with TruSight Oncology 500 application note. [illumina.com/content/dam/illumina-marketing/documents/products/appnotes/trusight-oncology-500-tmb-analysis-1170-2018-009.pdf](https://illumina.com/content/dam/illumina-marketing/documents/products/appnotes/trusight-oncology-500-tmb-analysis-1170-2018-009.pdf). Accessed September 20, 2023.

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© 2024 Illumina, Inc. All rights reserved.